

IRRADIATED RED BIOHAZARD BAG

Product Specification

Product Name:	Irradiated Red Biohazard Bag
Manufacturer Part Numbers:	BL2-REIRBB, BL2-3440REIRBB, BL2-25RENPIRBB
Description:	Irradiated red LLDPE biohazard bags with black printed Biohazard symbol on BL2-REIRBB and BL2-3440REIRBB. BL2-25RENPIRBB is a red LLDPE bag with no printing. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility Meets ASTM D1709 10 ⁻⁶ Sterility Assurance Level In-Compliance with AAMI Guidelines for Sterilization
Physical Data:	
Materials of Construction:	Red Linear Low-Density Polyethylene
Dimensions:	Thickness: 2 mil BL2-25REIRBB and BL2-25RENPIRBB: 25 in x 35 in (64 cm x 89 cm) BL2-3440REIRBB: 34 in x 40 in (86 cm x 102 cm)
For Use with: (as needed)	BLSSH series BinStand or similar container
Autoclavability:	Not Autoclavable
Chemical Compatibility:	N/A
Packaging:	5 Biohazard Bags per clear double bag pack, 10 packs per case
Traceability:	Lot numbers printed on each pack and each case Irradiation indicator on each pack and each case
Shelf Life:	36 months
Irradiation Standards:	Minimum irradiation dose 25kGy
Certificates:	Certification of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.