



Your strategic supplier for critical cleaning

CLEANROOM TAPES FOR A RANGE  
OF APPLICATIONS AND FUNCTIONS



*Excellence in cleaning since 1984*

Find it at [fishersci.com](http://fishersci.com) and [fishersci.ca](http://fishersci.ca)

 **fisher** scientific  
part of Thermo Fisher Scientific



**M**icronova is well equipped to meet your needs for cleanroom tape with a range of pressure sensitive and barrier tapes on a plastic core. The breadth of tapes includes low residue, low outgassing tapes for “clean” applications, heavy

duty tapes for clean build construction, high and low temperature tape, and ESD conductive grid, anti-static and polyimide tapes, as well as tapes for sealing gloves to garments.

## Vinyl Tapes



**Cleanroom Vinyl Tape**  
for the most critical cleanroom needs



**Irradiated Cleanroom Vinyl Tape**, individually double-bagged



**Economy Grade Cleanroom Vinyl Tape**, EGV for general controlled environment use

Micronova's cleanroom vinyl tapes are available in multiple widths and colors, ideal for production color marking. For the most critical applications, **Cleanroom Vinyl Tape** is best, while **Economy Grade Cleanroom Vinyl Tape** is used in a number of routine applications throughout controlled environments and **Irradiated Cleanroom Vinyl Tape** is designed for aseptic areas. Vinyl tapes are easily conformable, resistant to most common solvents and acids and hold up to autoclave temperatures when applied to a glass or stainless steel surface.

## Area Designation, Construction, Heavy Duty Tapes



**Cleanroom Construction Tape**  
Serrated Edge, Easy to Tear, High tack polyethylene construction and maintenance tape with rubber adhesive.



**MegaTape™**  
Medium density polyethylene tape, with clean peel medium tack acrylic adhesive. Low outgassing and excellent conformability




**Striped Safety Tape**  
Striped tape used for delineating and color coding specific equipment and areas. Poly-laminated vinyl film comes in three colors, Black/Yellow, Black/White, Red/White



**Cleanroom Barrier Tapes (Film)**  
High visibility red and yellow polyethylene film with bold black letters. Red “DO NOT ENTER” or Yellow “CAUTION” in both English and Spanish.

## Polyethylene Tapes

Micronova's polyethylene films are well suited for use in labs and controlled environments. Polyethylene film is low outgassing, has excellent chemical resistance and conformability, low in halogens and leachable chlorides. The acrylic adhesive allows for a clean release. **PCX** is clean-processed for your more critical environments, **PCX-IR** is irradiated for aseptic use and **EGP** is our economy grade polyethylene ideal for controlled and cleanroom support areas.



**Polyethylene Cleanroom Tape**



**Irradiated Polyethylene Cleanroom Tape**



**Economy Grade Polyethylene Tape**



### **CuffStrip™ Tape**

Lined polyethylene tape pre-cut to 18" lengths, ideal for sealing gloves to garments which eliminates skin exposure in critical environments.



### **MicroSeal™**

Perforated every 3" for easy tear, MicroSeal is used for sealing gloves to garments or wrapping plates used for environmental monitoring

## Electrostatic Discharge (ESD) Tapes



### **Polyimide Tape**

Polyimide film with silicone adhesive for high and low temp applications, -100F to 500F. 825A has acrylic adhesive for electronic industry.



### **NovaStat™ AntiStatic Tape**

Anti-static cellophane tape used for sealing and packaging in ESD sensitive environments.



### **Conductive Grid Tape**

Anti-static tape with conductive grids embedded between two anti-static layers helps pull static charge away from substrate.




## Specialty Tapes



### Autoclave Indicator Tape

Environmentally friendly lead-free autoclave indicating tape wound on plastic core. Stays on surface through steam sterilization. RoHS compliant.



### Double Coated Hi-Tack Tape


Double coated polyester film, with liner. Hi-tack adhesive works well when bonding dissimilar materials and hard to hold surfaces

## For Industry Uses See Micronova Application Flyers


(Click or Tap to View)



**Clinical & Research Laboratory**  
Cleaning and containment products for the Clinical Laboratory



**Cannabis Cleaning Products**  
Cleaning tools for specialized support in cannabis production areas



**Compounding Pharmacy**  
Cleaning and containment products for the Compounding Pharmacy



**High Technology Brochure**  
Products that support and simplify contamination control within high tech production



**Cleanroom Essentials**  
Truly essential products to keep your lab clean, safe, and efficient



**Lab, Production Getting Started**  
Specialty items that solve common challenges in labs and workspaces

Contact us today:

In the United States

Order online: [fishersci.com](http://fishersci.com)

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: [fishersci.ca](http://fishersci.ca)

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437

 **fisher scientific**  
part of Thermo Fisher Scientific