

Small Tools for Small Spaces



An efficient and ergonomic sanitization process is made simpler with the Micronova family of Small Tools for Small Spaces, designed for controlled environments and critical laboratory spaces. These tools are ideal for supporting contamination control in Isolators, hoods, RABs, Biosafety Cabinets, pass throughs, and work spaces.

All mop covers and small tools are cleanroom compatible, ISO 5 to 9.

- Covers are available in multiple materials from smooth polyester to textured foam
- Use with interchangeable fixed and extended handles
- Washable and reusable
- Holds up to strong cleaning solutions and disinfectants
- Covers available gamma irradiated.
- Certificates and Lot Traceability with each cover and tool
- Autoclavable

The Micronova Cleanroom **Slim T™ Mop** easily slips on to the 180-degree swivel T-Mop™ adapter for simple and thorough equipment cleaning. Stainless steel and acetal construction.



T-Mop - 10 inch

| | | |
|---|-----------|-----------|
| Stainless Steel T-Mop, 10 in, 25 cm | QDT-10 | 19035363D |
| Stainless Steel, 12 in, 30 cm Handle | SSU-9 | 191600118 |
| NovaPoly™ Slim T-Mop, 10 in, 25 cm | T6-10S | 19116693 |
| Polysorb™ Slim T-Mop, no foam, 10 in, 25 cm | T11-10SNS | 19166574 |

The Micronova **Hand Held Mop** conforms to the hand, then smoothly glides over surfaces for a complete clean or sanitization.

- Polyketone version is very durable, autoclavable
- Polypropylene version is reusable and then recyclable
- Fits 6" FlatMops



| | | |
|----------------------------------|-----------|----------|
| Polyketone Hand Held Mop | HHT-06PK | 19116722 |
| Polypropylene Hand Held Mop | HHT-06PP | 19116723 |
| NovaPoly™ Flatmop Cover, no foam | MZSM6-6NS | 19116724 |
| PolyMesh™ Flatmop Cover | MZSM68-6S | 19116725 |

The Micronova **Isolator Tool™ Mop** covers easily slide on to the unique adjustable 'swivel' joint tool allowing the ultimate in both ease and control of flat surface cleaning.



Isolator Tool & Covers - 9 inch

| | | |
|---|------------|----------|
| Plastic Isolator Tool, 9 in, 22 cm | QPIT-09 | 19435082 |
| Plastic Handle, 12 in, 30 cm | STAH-12WH | 18000400 |
| NovaPoly™ Irradiated Cover, 9 in, 22 cm | ITC6-9NSIR | 19166385 |
| Polysorb™ Isolator Tool Cover, no foam, 9 in, 22 cm | ITC11-9NS | 19435083 |

Micronova's **SlimLine™ Plastic Mop Adapter** is a polypropylene flat head mop frame adapter with 360° swivel joint and attachment site for handles that can be used with a range of Micronova flat head mop covers. Ideal when a standard mop head is too long for your cleaning application.

| | | |
|--|-------------|----------|
| Plastic SlimLine™ Mop, 9 in 22 cm | QPSL-09 | 19263611 |
| Polysorb™ Mop Cover, Irradiated 9 in, 22 cm | MZSM11-9SIR | 19263612 |
| NovaPoly™ Mop Cover, 9 in, 22 cm | MZSM56-9 | 19305730 |
| Extendible Plastic Handle, 18 to 36 in, 45-90 cm | STAEH-036WH | 18000411 |

A polypropylene diamond shaped flat head mop frame, the Micronova **Double Corner Tool**, with 360 swivel joint and attachment site for handles is perfect for cleaning coving, corners, pass throughs and small chambers.



| | | |
|---|-------------|----------|
| Plastic Double Corner Tool, 10 in, 25 cm | QPCCD-410 | 19116746 |
| NovaPoly™ Mop Cover, 10 in, 25 cm | CCDM6-410 | 19116753 |
| NovaPoly™ Mop Cover, Irradiated, 10 in, 25 cm | CCDM6-410IR | 19116792 |
| Plastic Threaded Handle, 24 in, 60 cm | STAH-24WH | 17020103 |

Contact us today:

In the United States

Order online: fishersci.com
Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Fax an order: 1-800-463-2996
Call customer service: 1-800-234-7437

