

## SPONGE, NOVACEL

### Product Specification

Product Name:	Sponge, NovaCel
Manufacturer Part Number(s):	SP4-533
Description:	A solid white absorbent PVA polyvinyl acetate sponge for cleaning small parts, equipment, and other lab applications. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction: Outer Layer	PVA polyvinyl acetate sponge
Dimensions: Width x Length x Height	SP4-533: 5 in x 3.3 in x 0.33 in (12.7 cm x 8.4 cm x .84 cm)
For Use with: (as needed)	N/A
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions
Packaging:	Individually bagged; 24 sponges per case.
Traceability:	Lot numbers printed on each bag and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website