

POLYMESH WIPE

Product Specification

Product Name:	PolyMesh Wipe
Manufacturer Part Number(s):	WP68-99H
Description:	A hemmed polyester mesh wipe with double ply polyester lining for abrasion cleaning and removal of tough soils. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer layer	100% Polyester braided mesh
Lining	100% knit polyester
Dimensions:	
Length x Width	9 in x 9 in (22.8 cm x 22.8 cm)
Thickness	1/16 in (1.5 mm)
For Use with: (as needed)	Abrasion cleaning to remove tough soils
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions
Packaging:	10 wipes per bag; 100 wipes per poly lined case
Traceability:	Lot numbers printed on each bag and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website