

STERILE BUCKETLINER

Product Specification

Product Name:	Sterile BucketLiner
Manufacturer Part Number(s):	BKL7-3119IR, BKL1-1818IR, BKL-2424IR, BKL-2430IR, BKL99-3521IR, BKL1-1818IRB, BKL993521IRB, BKL-2424IRBB, BKL-2424IRTBB
Description:	Sterile LDPE/Nylon laminate liners from FDA compliant LDPE with a horizontal orientation to be used as an impermeable barrier for lining buckets and containers. Can also be used to hold liquids going through the autoclave. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization.
Physical Data:	
Materials of Construction:	Clear Laminate: (LDPE) Polyethylene 75%/Cast Nylon 25%
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	
Thickness:	3 mil (0.076 mm)
Open Width x Height:	BKL7-3119IR: 31 in x 19 in (79 cm X 48 cm) BKL1-1818IR, -1818IRB: 18 in x 18 in (45 cm x 45 cm) BKL-2424IR, -2424IRBB, -2424IRTBB: 24 in x 24 in (60 cm X 60 cm) BKL-2430IR: 24 in x 30 in (60 cm x 90 cm) BKL99-3521IR, BKL993521IRB: 35 in x 21 in (89 cm x 53 cm)
For Use with: (as needed)	Buckets and containers similar in shape to: BKL7-3119IR: B-7 buckets BKL1-1818IR, -1818IRB: 1 gal pails BKL-2424IR, -2424IRBB, -2424IRTBB: 6-gal round or square buckets BKL-2430IR: 8-gal round bucket BKL993521IRB, -3521IR: B-99 bucket
Autoclavability:	Autoclavable – Will hold up to autoclaving with liquid contents
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions; avoid chlorites

Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	BKL7-3119IR, BKL-2424IR, BKL-2430IR, BKL99-3521IR: Individually double bagged; 24 liners per case BKL1-1818IR: Individually double bagged; 50 liners per case BKL1-1818IRB: 10 liners per pack; double bagged; 5 packs per case. 50 liners total BKL993521IRB: 6 liners per pack; double bagged; 4 packs per case. 24 liners total BKL-2424IRBB: 6 liners per pack, double bagged; 8 packs per case. 48 liners total BKL-2424IRTBB: 1 pack of 6 liners; single bagged; 4 packs of 6 liners; double bagged; 2 packs of 24 liners; triple bagged. 48 liners total
Traceability:	Every pack and case is labeled with its lot number and part number
Shelf Life:	36 months
Irradiation Standards:	Minimum irradiation dose 25kGy
Certificates:	Certificate of Irradiation and Certificate of Conformance shipped with each order and available on website.