

STERILE BUCKETLINER

Product Specification

Product Name:	Sterile BucketLiner
Manufacturer Part Numbers:	BKL7-3119IR, BKL1-1818IR, BKL-2424IR, BKL-2430IR, BKL99-3521IR, BKL993521IRB, BKL-2424IRBB, BKL-2424IRTBB
Description:	LDPE/Nylon blend liner for buckets and storage containers. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization.
Physical Data:	
Materials of Construction:	Clear LDPE 75%, nylon 25% laminate
Dimensions:	3 mil
Thickness:	BKL7-3119IR: 31 in x 19 in (79 cm X 48 cm)
Open Width x Height:	BKL1-1818IR: 18 in x 18 in (45 cm x 45 cm)
	BKL-2424IR, BKL-2424IRBB, BKL-2424IRTBB: 24 in x 24 in (60 cm X 60 cm)
	BKL-2430IR: 24 in x 30 in (60 cm x 90 cm)
	BKL99-3521IR, BKL993521IRB: 35 in x 21 in (89 cm x 53 cm)
For Use with: (as needed)	Buckets and containers similar in shape to: BKL7-3119IR: B-7 buckets BKL1-1818IR: 1 gal pails BKL-2424IR, BKL-2424IRBB and BKL-2424IRTBB: C-2PW buckets BKL-2430IR: 8-gal round bucket BKL99-3521IR and BKL993521IRB: B-99 bucket
Autoclavability:	Autoclavable - Will hold up to autoclaving with liquid contents
Chemical Compatibility:	Avoid Chlorites
Packaging:	BKL7-3119IR, BKL-2424IR, BKL-2430IR, BKL99-3521IR: Individually double bagged, 24 bags per case BKL1-1818IR: Individually double bagged, 50 bags per case BKL993521IRB: 6 Liners per Pack, Double Bagged, with 4 pkgs per case. 24 bags per case total BKL-2424IRBB, BKL-2424IRTBB: 6 Liners per Pack, Double Bagged, with 8 pkgs per case. 48 bags per case total
Traceability:	Lot numbers printed on each package and each case
Shelf Life:	36 Months
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation and Certificate of Conformance shipped with each order and available on website.