

STERILE BUCKETLINER™ WITH EXPIRATION DATE

Product Specification

Product Name:	Sterile BucketLiner with Expiration Date
Manufacturer Part Number(s):	BKL7-3119IRX
Description:	Sterile, LDPE/Nylon laminate liners with a labeled expiration date from FDA compliant LDPE with a horizontal orientation to be used as an impermeable barrier for lining buckets and containers. Can also be used to hold liquids going through the autoclave. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization
Physical Data:	
Materials of Construction:	Clear Laminate: (LDPE) Polyethylene 75%/Cast Nylon 25%
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Thickness: 3 mil (0.076 mm)
	Open Width x Height: 31 in x 19 in (79 cm x 48 cm)
For Use with: (as needed)	Buckets and containers similar in shape to: BKL7-3119IRX: B-7/B-7AC buckets; BucketBinder™ Multi Bucket Systems
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with standard disinfectants and cleaning solutions; avoid chlorites
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources
Packaging Configuration:	Individually double bagged; 24 bags per case; using a ½-inch white tape seal on outer poly bag liner.
Traceability:	Every pack and case is labeled with its lot number and part number



Revision 001

Shelf Life:	36 months
Irradiation Standards:	Minimum irradiation dose 25kGy
Certificates:	Certificate of Irradiation and Certificate of Conformance shipped with each order and available on website.