

# BUCKETLINER

## Product Specification

Product Name:	BucketLiner
Manufacturer Part Number(s):	BKL7-3119, BKL-2424, BKL1-1818, BKL99-3521
Description:	LDPE/Nylon laminate liners from FDA compliant LDPE with a horizontal orientation to be used as an impermeable barrier for lining buckets and containers. Can also be used to hold liquids going through the autoclave. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Clear Laminate: (LDPE) Polyethylene 75%/Cast Nylon 25%
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	
Thickness:	3mil (0.076 mm)
Open Width x Height:	BKL7-3119: 31 in x 19 in (79 cm x 28 cm) BKL-2424: 24 in x 24 in (60 cm x 60 cm) BKL1-1818: 18 in x 18 in (45 cm x 45 cm) BKL99-3521: 35 in x 21 in (89 cm x 53 cm)
For Use with: (as needed)	Buckets and containers similar in shape to: BKL7-3119: B-7/B-7AC series buckets; BucketBinder™ Multi Bucket Systems BKL-2424: 6-gal round or square buckets; C-2PW BKL1-1818: 1-gal pails BKL99-3521: B-99 stainless steel bucket
Autoclavability:	Autoclavable – Will hold up to autoclaving with liquid contents
Chemical Compatibility:	Good compatibility with standard disinfectants and cleaning solution; avoid chlorites
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.

Packaging Configuration:	BKL7-3119, BKL-2424 and BKL99-3521: 20 liners per pack; 10 packs per case; 200 liners total BKL1-1818: 50 liners per pack; 3 packs per case; 150 liners total
Traceability:	Every pack and case is labeled with its lot number and part number
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available on request or website