

## STERILE BINLINER, 2 MIL

### Product Specification

Product Name:	Sterile BinLiner, 2 mil
Manufacturer Part Number(s):	BL2-3219IRBB, BL2-3219PCIR, BL2-3440IRTBB, BL2-3625IRBB, BL2-40TRIR, BL2-40TRIRB, BL2-40TRIRBB, BL2-40TRPCIR, BL2-1016IRBB
Description:	Irradiated polyethylene liners from FDA compliant LDPE. May be used for material disposal or transfer or as a lightweight protective equipment cover. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level In-compliance with AAMI Guidelines for Sterilization.
Physical Data:	
Materials of Construction:	Clear Low-Density Polyethylene, Meets FDA Regulation 21 CFR, 177.1520, 3.1
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of 100% recycled material
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal and human origin.
Dimensions:	Thickness: 2 mil Open Width x Height: BL2-3219 series: 32 in x 19 in (81 cm x 48 cm) BL2-3440IRTBB: 34 in x 40 in (86.4 cm x 101.6 cm) BL2-3625IRBB: 36 in x 25 in (90 cm x 63 cm) BL2-40TR series: 40 in x 48 in (100 cm x 120 cm) Gusset: 23 in x 17 in (58 cm x 43 cm) BL2-1016IRBB: 10 in x 16 in (25.4 cm x 40.6 cm)
For Use with: (as needed)	N/A
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common cleaners and solutions.
Storage Conditions:	Ambient temperature (Recommended between 15°C - 25°C). Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet light sources.
Packaging Configuration:	B Suffix: 1 liner per bag, 5 liners per pack, 5 liners per double bag BB Suffix: BL2-1016IRBB: 10 liners per double bag TBB Suffix: BL2-3440IRTBB: 50 liners per case, triple bagged with additional case liner. PC Suffix: Roll of liners per double bag
Traceability:	Lot numbers printed on each package and each case



Revision 003

Shelf Life:	36 months
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.