

# RED BIOHAZARD BAG, 2 MIL

## Product Specification

Product Name:	Red Biohazard Bag, 2 mil
Manufacturer Part Number(s):	BL2-3440REBB
Description:	<p>A red biohazard polyethylene liner made from FDA compliant LDPE with "Biohazard" and a Biohazard symbol printed on each bag. Biohazard bags are designed for the collection and containment of biomedical waste.</p> <p>Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.</p>
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Red Low-Density Polyethylene. Meets FDA Regulation 21 CFR, 177.1520, 3.1
Dimensions:	34 in x 40 in (86 cm x 100 cm)
For Use with: (as needed)	BLSSH-1515 or similar container
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common cleaners and solutions.
Packaging:	5 liners per double bag
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Not irradiated
Certificates:	Certificate of Conformance available on request and on website.