

## SMALL BINLINER W/ ZIP TIES, IRRADIATED, 2MIL

### Product Specification

Product Name:	Small BinLiner w/ zip ties, Irradiated, 2mil
Manufacturer Part Number(s):	BL2-1016ZIRBB
Description:	Irradiated polyethylene liners from FDA compliant LDPE with irradiated nylon zip ties. May be used for material disposal or transfer. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level In-compliance with AAMI Guidelines for Sterilization
Physical Data:	
Materials of Construction:	
Liners:	Clear Low-Density Polyethylene, Meets FDA Regulation 21 CFR, 177.1520, 3.1
Zip Ties:	Nylon
Dimensions:	
Thickness:	2 mil
Open Width x Height:	10" x 16" (25 cm x 40.6 cm)
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common cleaners and solutions
Packaging:	
BB Suffix:	10 bags and zip ties per pack double-bagged, 20 packs to a case. 200 total bags
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.