

BINLINER, 4 MIL

Product Specification

Product Name:	BinLiner, 4 mil
Manufacturer Part Number(s):	BL4-1620B, BL4-3440B, BL4-3440BB, BL4-4060BB, and BL4-40TRBB
Description:	Polyethylene liners from FDA compliant LDPE. May be used for material disposal or transfer or as a lightweight protective equipment cover. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Clear Low- Density Polyethylene Meets FDA Regulations 21CFR177.1520, 3.1
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of 100% recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origins.
Dimensions:	Thickness: 4 mil BL4-1620B: 16 in x 20 in (40 cm x 50 cm) No Gusset BL4-3440B, -3440BB: 34 in x 40 in (85 cm x 100 cm) Gusset Dimensions: 20 in x 14 in (20 cm x 35 cm) BL4-4060BB: 40 in x 60 in (100 cm x 150 cm) Gusset Dimensions: 23 in x 17 in (58 cm x 43 cm) BL4-40TRBB: 40 in x 48 in (100 cm x 120 cm) Gusset Dimensions: 23 in x 17 in (58 cm x 43 cm)
For Use with: (as needed)	BinLiner Stands: BL4-3440 series use BLSSH-1515 BL4-40TR and BL4-4060 series use BLSSH-1624
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	B Suffix: 1 liner per bag, 50 per case BB Suffix: 5 liners per double bag, 50per case
Traceability:	Lot numbers and part numbers printed on each pack and each case



Revision 000

Shelf Life:

N/A

Irradiation Standards:

Not Irradiated

Certificates:

Certificate of Conformance is available on request and available on website.