

## STERILE BINLINER, 4 MIL

### Product Specification

Product Name:	Sterile BinLiner, 4 mil
Manufacturer Part Number(s):	BL4-1620IRBB, BL4-3440IRBB, BL4-3440PCIR, BL4-3440IRTB, BL4-4060IRBB, BL4-40TRIR, BL4-40TRIRB, BL4-40TRIRBB
Description:	Irradiated polyethylene liners from FDA compliant LDPE. May be used for material disposal or transfer or as a lightweight protective equipment cover. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10 <sup>-6</sup> Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization.
Physical Data:	
Materials of Construction:	Clear Low- Density Polyethylene, Meets FDA Regulation 21CFR177.1520, 3.1
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Thickness: 4 mil Open Width x Height: BL4-1620IRBB: 16 in x 20 in (40 cm x 50 cm) Gusset Dimensions: No gusset Open Width x Height: BL4-3440IRBB, -3440PCIR, -3440IRTB: 34 in x 40 in (85 cm x 100 cm) Gusset Dimensions: 20 in x 14 in (50 cm x 35 cm) Open Width x Height: BL4-4060IRBB: 40 in x 60 in (100 cm x 120 cm) Gusset Dimensions: 23 in x 17 in (58 cm x 43 cm) Open Width x Height: BL4-40TRIR, -40TRIRB, -40TRIRBB: 40 in x 48 in (100 cm x 120 cm) Gusset Dimensions: 23 in x 17 in (58 cm x 43 cm)
For Use with: (as needed)	BinLiner Stands: BL4-3440 series use BLSSH-1515 BL4-40TR and BL4-4060 series use BLSSH-1624
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions.
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.

Packaging Configuration:

B Suffix: 1 liner per bag, 5 bags per pack  
BB Suffix: 5 liners per double bag; with the exception of BL4-1620IRBB which is 10 liners per double bag (smaller bag)  
PC Suffix: 50 liners on a roll double bagged.  
TB Suffix: 5 liners per bag, triple bagged, 10 pkgs per two additional bags sealed with cleanroom tape.

Traceability:

Lot numbers and part numbers are printed on each pack and each case.

Shelf Life:

36 months

Irradiation Standards:

Minimum irradiation dose 25kGy

Certificates:

Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.