

STERILE YELLOW BIOHAZARD BAG, 4 MIL

Product Specification

Product Name:	Sterile Yellow Biohazard Bag, 4 mil
Manufacturer Part Number(s):	BL4-32YEIRBB
Description:	An irradiated yellow biohazard polyethylene liner made from FDA compliant LDPE with "Biohazard" and a Biohazard symbol printed on each bag. Biohazard bags are designed for the collection and containment of biomedical waste. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility 10-6 Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization
Physical Data:	
Materials of Construction:	Yellow Low-Density Polyethylene with black printed BioHazard symbol, Meets FDA Regulation 21 CFR, 177.1520, 3.1
Dimensions:	
Thickness:	4 mil
Open Width X Height:	32 in X 40 in (80 cm x 100 cm)
For Use with: (as needed)	BLSSH-1515 BinLiner Stand
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common cleaners and solutions.
Shelf Life	36 Months
Packaging:	5 liners per double bag
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website. Validation paperwork available upon request.