

STERILE BINLINER, 6 MIL

Product Specification

Product Name:	Sterile BinLiner, 6 mil
Manufacturer Part Number(s):	BL6-4048IRBB; BL6-4048PCIR
Description:	Irradiated polyethylene liners from FDA compliant LDPE. May be used for material disposal or transfer or as a lightweight protective equipment cover. Manufactured at Micronova's Torrance, CA facility Designed for clean and controlled environments, ISO5 to ISO9
Manufacturing Standards:	ISO 9001 Certified Facility 10 ⁻⁶ Sterility Assurance Level per AAMI TIR33:2005
Physical Data:	
Materials of Construction:	Clear Low-Density Polyethylene, Meets FDA Regulation 21 CFR, 177.1520, 3.1
Dimensions:	
Thickness:	6 mil
Open Width x Height	40 in x 48 in (100 cm x 120 cm);
Gusset Dimensions	no gusset
For Use with: (as needed)	BLSSH-1624 BinLiner Stand C-2BAG cart for rolls
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common cleaners and solutions
Shelf life:	36 months
Packaging:	
BB Suffix	5 liners per double bag
PC Suffix	Roll of liners per double bag
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Minimum irradiation dose 25kGy
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.