

# BINLINER STAND

## Product Specification

Product Name:	BinLiner Stand
Manufacturer Part Number(s):	BLSSH-1624, BLSSH-168, BLSSH-1515, BLSSH-1810, BLSSH-46
Description:	Stainless steel frame for use with Micronova's series of BinLiners. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	304 Stainless Steel 1/4 in Rods – Electropolished Polypropylene Connectors
Packaging Materials:	Packaged in LDPE bag
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Width x Length x Height: BLSSH-1624: 16 in x 25 in x 34 in (40.6 cm x 63.5 cm x 86.4 cm) BLSSH-168: 16 in x 8 in x 28 in (40.6 cm x 20.3 cm x 71.1 cm) BLSSH-1515: 15 in x 15 in x 28 in (38.1 cm x 38.1 cm x 71.1 cm) BLSSH-1810: 18 in x 10 in x 30 in (45.7 cm x 25.4 cm x 76.2cm) BLSSH-46: 5 in x 5 in x 11 in (12.7 cm x 12.7 cm x 27.9 cm)
For Use with: (as needed)	BinLiner Bags
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (Recommended between 15°C - 25°C). Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet light sources.
Packaging Configuration:	Individually bagged
Traceability:	Lot numbers and part numbers printed on each pack and each case
Shelf Life:	N/A
Irradiation Standards:	N/A
Certificates:	Certificate of Conformance available upon request