

PLASTIC CABLE TIES, IRRADIATED

Product Specification

Product Name:	Plastic Cable Ties, Irradiated
Manufacturer Part Numbers:	BLT-12IR, BLT-12TIR
Description:	12" Long solid nylon tie, irradiated. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Nylon
Dimensions:	
Width:	23 in (.58 cm)
Length:	11.75 in (30 cm)
Thickness:	.05 in (.13 cm)
For Use with: (as needed)	All Micronova Bin and Bucket Liners
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	BLT-12IR: 5 each per Bag, 5 packs per Double Bag, 4 packs of 25 per case. Total 100/cs BLT-12TIR: 5 each per bag, 25 packs of 5 triple bagged. 5 packs of 25 per case. Total 125/cs.
Traceability:	Lot number and part number printed on each bag and each case
Shelf Life:	N/A
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificates of Conformance and Certificate of Irradiation Processing shipped with each order and available on request.