

# AUTOCLAVE INDICATOR TAPE, LEAD FREE

## Product Specification

Product Name:	Autoclave Indicator Tape, Lead Free
Manufacturer Part Number(s):	AUTOPC-1, AUTOPC-12, AUTOPC-15, AUTOPC-2, AUTOPC-25, AUTOPC-3, AUTOPC-24
Description:	A lead-free indicator tape for steam sterilization cycles. Green indicator lines change to black at completion of the sterilization cycle in both gravity discharge and vacuum autoclave. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility RoHS compliant
Physical Data:	
Materials of Construction:	Film: Saturated crepe paper Adhesive: Synthetic elastomer Core: Plastic Color: Natural
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Core: 3 in (7.6 cm) Length: 60 yards (55 m) Width: AUTOPC-34: in (1.9 cm) AUTOPC-12: in (1.3 cm) AUTOPC-1: 1 in (2.5 cm) AUTOPC-15: 1.5 in (3.8 cm) AUTOPC-2: 2 in (5 cm) AUTOPC-25: 2.5 in (6.4 cm) AUTOPC-3: 3 in (7.6 cm)
For Use with: (as needed)	P-52, P-52W Tape Dispensers
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources. Store product on a flat/level surface to prevent telescoping of tape.



Revision 002

Packaging Configuration:	Individually bagged AUTOPC-1: 36 rolls per case AUTOPC-12: 16 rolls per case AUTOPC-15, AUTOPC-2: 24 rolls per case AUTOPC-25, AUTOPC-3: 12 rolls per case AUTOPC-34: 16 rolls per case
Traceability:	Part number and lot number printed on each bag and case
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance is available upon request and on website.