

ISOLATOR TOOL™ ADAPTER, PLASTIC

Product Specification

Product Name:	Isolator Tool™ Adapter, Plastic
Manufacturer Part Number(s):	QPIT-04, QPIT-09
Description:	Plastic bladed Isolator Tool with a flat, angled head and a 360° swivel joint lockable at 180° for optimized cleaning. Features a threaded attachment compatible with Micronova® STA series plastic handles. Designed for controlled and cleanroom environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Plate: Acetal Plastic Swivel: Polyketone Handle Connection: Polyketone
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Height: Swivel flat: 2 ½ in (6.4 cm) Swivel extended: 5 ½ in (14 cm) Plate: QPIT-04: 2 x 4 in (5.1 x 10.2 cm) QPIT-09: 2 x 9 in (5.1 x 22.9 cm) Weight or Capacity: Less than 1 pound with cover included
For Use with: (as needed)	Isolator Tool™ covers and threaded STA Series fixed length or extendable handles
Autoclavability:	Autoclavable (Recommended up to 5 cycles). Before autoclaving, rinse thoroughly to remove disinfectant residue. Loosen collar on adapter stem prior to autoclaving. Autoclaving unit with collar locked down can cause threads to freeze in place. Store flat and avoid manipulation immediately after autoclave.
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (Recommended between 15°C - 30°C.)



Revision 002

Packaging Configuration:

Individually bagged

Traceability:

Part Number and Lot Number on each bag

Shelf Life:

N/A

Irradiation Standards:

Not Irradiated

Certificates:

Certificate of Conformance available upon request and/or website.