

MZONE FLATMOP COVER, NOVAPOLY, NO FOAM

Product Specification

Product Name:	MZone FlatMop Cover, NovaPoly, No Foam
Manufacturer Part Number(s):	MZSM6-6NS, MZSM6-9NS, MZSM6-9NSBB
Description:	A slip-on polyester mop cover, with no foam for cleaning small equipment. Designed for isolators, bio safety cabinets, and other small clean and controlled environments. Manufactured at Micronova's Torrance, CA facility.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	100% Polyester, 2 ply
Dimensions:	MZSM6-6NS: 4 in X 6 in (10.2 cm x 15.2 cm) MZSM6-9NS, MZSM6-9NSBB: 4 in X 9 in (10.2 cm x 22.9 cm)
For Use with: (as needed)	MZSM6-6NS: Use with Handheld Tools, either HHT- 06PK or HHT-06PP MZSM6-9NS: Use with QDSL Stainless Steel SlimLine mop adapter, SSU handles or QPSL Plastic SlimLine mop adapter, STA series handles
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Shelf Life:	NA
Packaging:	Individually bagged. MZSM6-6NS: 6 bags/pack, 4 packs/case MZSM6-9NS: 24 mops/case MZSM6-9NSBB: 10ea per pack, 20 packs per case, 200ea per case total.
Traceability:	Lot numbers and part numbers printed on each pack and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available on request