

HANDHELD SCRAPER TOOL

Product Specification

Product Name:	Handheld Scraper Tool
Manufacturer Part Number(s):	HHSWH-8 and HHSWH-10
Description:	White handheld scraper tool for removing residues and build up. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	One piece, injection molded polypropylene construction. Made from FDA approved materials. No seams or cracks for bacterial growth. Will not scratch metal surfaces.
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of 100% recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal and human origin.
Dimensions:	HHSWH-8: 3 in X 8 in (8cm x 20 cm) white scraper HHSWH-10: 4.4 in X 9.7 in (11 cm x 25 cm) white scraper
For Use with: (as needed)	N/A
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with common cleaning solutions and disinfectants.
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	Each scraper is single bagged
Traceability:	Lot numbers and part numbers printed on each bag and each case
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request