

AQUAHOL STERILE ALCOHOL

Product Specification

Product Name:	AquaHol Sterile Alcohol
Manufacturer Part Number(s):	IC-220
Description:	Sterile hand and glove sanitizer comprised of 70% isopropyl alcohol with 30% WFI quality purified water, filtered to 0.1 microns and meets USP requirements. Ideal for use over gloves prior to entering into aseptic areas and between critical processes. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 thru 9.
Manufacturing Standards:	ISO 9001 Certified Facility. 10-6 Sterility Assurance Level per VDMax Sterilization Method.
Physical Data:	
Materials of Construction:	Sanitizer: 70% Isopropyl alcohol with 30% purified water/WFI quality, filtered to 0.1 microns. Container: HDPE
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Net contents: 500mL
For Use with: (as needed)	For use with the ICDISP-3 Touch-free MicroDispenser
Autoclavability:	Not Autoclavable
Chemical Compatibility:	The product is chemically stable. Under normal conditions of storage and use, hazardous re-actions will not occur.
Storage Conditions:	Ambient temperature (recommended between 15°C - 25°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources. Keep away from heat/sparks/open flames/hot surfaces – No smoking.
Packaging Configuration:	6 bottles per poly-lined case.
Traceability:	Lot numbers printed on each case and on each individual bottle.
Shelf Life:	2 years after manufacture date
Irradiation Standards:	Irradiated



Revision 002

Certificates:

Certificate of Analysis, Certificate of Sterility and Certificate of Irradiation provided with each order and available upon request