

AQUAHOL STERILE ALCOHOL

Product Specification

Product Name:	AquaHol Sterile Alcohol
Manufacturer Part Number(s):	IC-220
Description:	Sterile hand and glove sanitizer comprised of 70% isopropyl alcohol with 30% WFI quality purified water, filtered to 0.1 microns and meets USP requirements. Ideal for use over gloves prior to entering into aseptic areas and between critical processes. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility. 10⁻⁶ Sterility Assurance Level per VDMax Sterilization Method
Physical Data:	
Materials of Construction:	
Sanitizer:	70% Isopropyl alcohol with 30% purified water/WFI quality, filtered to 0.1 microns.
Container:	HDPE
Dimensions:	
Net Contents	500 mL
For Use with: (as needed)	For use with the ICDISP-3 Touch-free MicroDispenser
Autoclavability:	Not Autoclavable
Chemical Compatibility:	The product is chemically stable. Under normal conditions of storage and use, hazardous reactions will not occur.
Shelf Life:	3 years after manufacture date
Packaging:	6 bottles per poly-lined case
Traceability:	Lot numbers printed on each case and on each individual bottle
Irradiation Standards:	Irradiated
Certificates:	Certificate of Analysis provided with each order and available upon request