

Micronova[®]

MANUFACTURING INC.

Your strategic supplier for critical cleaning

PENS AND MARKERS FOR
SUPERIOR DOCUMENTATION



Excellence in cleaning since 1984

Meeting requirements for Good Documentation Practices is imperative across controlled environments. From low sodium ink pens for Electronics/Semiconductor paper records to irradiated markers for labeling plastic, glass and paper in aseptic areas and RABS, we have the pens and packaging to meet your requirement.



Cleanroom Pens

Cleanroom Pens feature a low-sodium, indelible ink and a special high-density poly barrel to limit particle generation.

Available in Black, Blue and Red.

Certified for Class 100.



Cleanroom Click Pens

Cleanroom Click Pens have the same clean, low sodium, indelible ink features as PEN-01 and PEN-02 (above) with an easy "click" retractable nib.

Available in black and blue ink.

Cleanroom Pens, Cleanroom Click Pens, IPA Resistant LabMarkers and Rollerball GelPens are available gamma irradiated – suited for Aseptic ISO 5, and other clean environments.



Rollerball IPA Resistant GelPen – Irradiated

A fine point, archival ink gel pen for paper records.

Available in blue ink.



IPA-Resistant LabMarkers™

Lab Markers feature permanent, quick drying ink that leaves a clear impression on plastic bags or disposable labware. Has resistance to IPA overspray and smearing. Can be washed off with soap and water.



Sharpie® Fine Point Lab Markers

The industry standard, pen style permanent marker. Available in original "capped" (above) and 'retractable' styles (below).



Sharpie® UltraFine Point Lab Markers

With the same industry standard permanent ink, this marker has an ultra fine tip for precise labeling and documentation. Capped. Available in blue or black ink.



UltraFine Smear-Proof Markers

A permanent marker that gives the thinnest permanent line at 0.3 mm. The quick dry ink works well with lab samples and slides, glass, plastic as well as paper. Long nib allows for writing in limited spaces.

Sharpie® Fine Point and Ultra Fine Point Markers are gamma irradiated – suited for Aseptic, ISO 5, and other clean environments.

Part No.

Low Sodium Ink

Rollerball Gel/ Archival Ink

Fiber Tipped "Marker" Pen

Capped Nib

Click/ Retractable

Blue Ink

Black Ink

Red Ink

Irradiated

IPA Resistant Ink

Fine Point

Ultra Fine Point

Pens Per First Bag

Pens Per Double Bag

Pens Per UOM Sold

Cleanroom Pens

PEN-01	✓			✓			✓						10	10	10
PEN-02	✓			✓		✓							10	10	10
PEN-03	✓			✓				✓					10	10	10
PEN-01IR	✓			✓			✓		✓				10	10	100
PEN-02IR	✓			✓		✓			✓				10	10	100
PEN-03IR	✓			✓				✓	✓				10	10	100

Cleanroom Click Pens

PEN-11	✓				✓		✓						10	10	10
PEN-12	✓				✓	✓							10	10	10
PEN-11IR	✓				✓		✓		✓				10	10	100
PEN-12IR	✓				✓	✓			✓				10	10	100

Rollerball IPA-Resistant GelPen - Irradiated

PEN-18IR		✓			✓	✓			✓	✓	✓		1	4	12
----------	--	---	--	--	---	---	--	--	---	---	---	--	---	---	----

IPA-Resistant LabMarkers™

PEN-40			✓	✓			✓			✓	✓		5	10	10
PEN-40IR			✓	✓			✓		✓	✓	✓		1	5	10

Sharpie® Fine Point Markers

PEN-20IR			✓	✓			✓		✓		✓		1	4	12
PEN-22IR			✓	✓		✓			✓		✓		1	4	12
PEN-23IR			✓	✓				✓		✓			1	4	12
PEN-25IR			✓		✓		✓		✓		✓		1	4	12
PEN-26IR			✓		✓	✓			✓		✓		1	4	12
PEN-27IR			✓		✓			✓		✓			1	4	12

Ultra Fine Point Markers

PEN-51			✓	✓			✓				✓		4	12	12
PEN-52			✓	✓		✓					✓		4	12	12
PEN-51IR			✓	✓			✓		✓		✓		1	4	12
PEN-52IR			✓	✓		✓			✓		✓		1	4	12
PEN-53IR			✓	✓				✓		✓			1	4	12
PEN-61IR			✓	✓		✓			✓		✓		1	4	12



Micronova Manufacturing, Inc.
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com

