

CLEANROOM CLICK PEN - IRRADIATED

Product Specification

Product Name:	Cleanroom Click Pen - Irradiated
Manufacturer Part Number(s):	PEN-11IR, PEN-12IR
Description:	A low-sodium ink, click ball point pen, designed with a special poly barrel to limit particle generation. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	White high-density poly barrel Low-sodium ink PEN-11IR: black ink PEN-12IR: blue ink
Dimensions:	6 in x 3/8 in
For Use with: (as needed)	Cleanroom documents
Autoclavability:	Not autoclavable
Chemical Compatibility:	Holds up to wipe down with IPA, bleach and phenols. Integrity of printed barrel will be compromised with thinners and strong solvents.
Shelf Life:	24 months
Packaging:	Packaged 10 per pack / 10 packs per case
Traceability:	Lot numbers printed on each outer bag and each case
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation and Certificate of Conformance shipped with each order and available on website.