

# Double CornerCleaner™



**M**icronova's **Double CornerCleaner™** is designed to easily clean a wide range of flat surfaces in the cleanroom. The diamond-shaped adapter is ideal for cleaning equipment corners, pass-throughs, BSCs, hoods, chambers and other small spaces.

The Double CornerCleaner Mop Covers are available with or without a foam interior in a variety of cleanroom fabrics to enhance the distribution of sanitizing agents and remove contaminants.

The covers are available gamma irradiated for aseptic applications. The tool's durable polypropylene base features a unique 360° swivel joint that can be locked at 180° for more controlled cleaning. The system is lightweight and easy to use for cleaning or disinfecting in ISO 5 to 9 areas.



## SPECIFICATIONS

- QD frame adapter for use with SSU series handles
- QP frame adapter for use with STA series ACME threaded handles
- 360° plastic swivel joints can lock at 180° for more control
- Adapters and covers are autoclavable
- Covers are available in a variety of fabrics for specialized cleaning
- Covers available gamma irradiated

### Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility
- Autoclavable
- Good chemical resistance to disinfectants
- Minimum Irradiation Dose: 25 kGy
- 10<sup>-6</sup> Sterility Assurance Level
- In compliance with AAMI Guidelines for Sterilization

### Materials of Construction:

QD Plate:	Polypropylene
QD Swivel Joint:	Polypropylene, Stainless Steel, Polyketone
QP Plate:	Polypropylene
QP Swivel Joint:	Polypropylene, Polyketone

### Covers:

Novapoly™	100% Smooth polyester knit provides excellent disinfectant dispersal
Microfiber	80% Polyester, 20% polyamide mop is highly absorbent and effective for wet and dry mopping

### Dimensions:

QD and QP Plate:	4 in x 10 in (10 cm x 25.4 cm)
Swivel Flat:	2-1/2 in (6 cm)
QD Swivel Extended:	9 in (23 cm)
QP Swivel Extended:	5-1/2 in (14 cm)

**Weight or Capacity:** 1 lb. (.45 kg)



**Micronova Manufacturing, Inc.**  
3431 W. Lomita Boulevard  
Torrance, CA 90505-5010 USA

Tel: 310/784-6990  
Fax: 310/784-6980

SalesInfo@micronova-mfg.com  
www.micronova-mfg.com