

T-Mop™ Adapter



The Micronova **T-Mop™** series of mop adapters feature a simple T shape with 180° pivot for angling the mop head. This allows the round stainless steel or square acetyl bar to simply roll down the wall, across the ceiling, or over the floor.

T-Mop adapters support T-Mop™ or Slim T™ mop covers, which easily snap onto the cross bar. The acetyl bar style works best with thick foam T-Mop or Slim T mop covers. Both the round stainless and square acetyl T-Mop adapters work with stainless steel quick disconnect handles (QD version) or with plastic threaded handles (QP version).



SPECIFICATIONS

Use with Micronova mops and hardware

- 180° pivot allows mop to clean easily under equipment or in awkward corners and tight spaces
- T-Mops work well on ceilings, walls, floors, soft walls and equipment exteriors
- Stainless-steel tube frame fits all T-Mop and Slim T covers
- QD, Quick Disconnect attachment to SSU series stainless handles
- QP, Threaded end for attachment to STA series plastic handles.
- Adapters are autoclavable and have good chemical resistance

Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility
- Autoclavable
- Good chemical resistance to disinfectants and cleaning solutions

Materials of Construction:

Mop Bar:

Standard	Stainless Steel Tube
Square	Acetal Plastic

180° Angle Joint:

Acetal Plastic

Handle Connector:

QD Series	Stainless Steel Acetal
QP Series	Plastic

Dimensions:

Mop Bar:

Length	10, 14, or 18 in (25, 35 or 45 cm)
--------	------------------------------------

Diameter:

Standard	1 inch (2.5 cm)
Square	0.75 in (2.0 cm)

Handle Connection:

QD Extension	4.5 in (11.43 cm)
QP Extension	2.5 in (6.35 cm)

Weight or Capacity: Approximately 1 pound



Micronova Manufacturing, Inc.
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com