

IRRADIATED BAG CUTTER

Product Specification

Product Name:	Irradiated Bag Cutter
Manufacturer Part Number(s):	BC-01IR and BC-02IR
Description:	Bag cutters are available in two styles. The orange style bag cutter offers an automatically retractable safety blade guard with fold back cover and the yellow hook style. Cutters are disposable with no blade changes required for either style. Excellent for opening raw material packages. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	BC-01IR: Orange, plastic with retractable blade BC-02IR: Yellow, plastic with hook style safety cutter
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of 100% recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origins.
Dimensions:	BC-01IR: Length: 6 in (15.2 cm) Blade head width: 1 ½ in (3.8 cm) Handle width: 1 1/8 in (2.9 cm) Height: ½ in (1.8 cm) BC-02IR: Length: 6 ½ in (16.5 cm) Blade head width: 2 ½ in (6.3 cm) Handle width: 1 1/8 in (2.9 cm) Height: 7/16 in (1.1 cm)
For Use with: (as needed)	N/A
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C). Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.



Revision 001

Packaging Configuration:	2 each double bagged pack, 25 packages per case, 50 bag cutters total.
Traceability:	Part number and Lot number printed on each bag and each case
Shelf Life:	60 Months
Irradiation Standards:	Irradiated
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.