

SealTech™ Resealable Device Pouch, Irradiated

Product Specification

Product Name:	SealTech™ Resealable Device Pouch, Irradiated
Manufacturer Part Number(s):	BKPH-58IRB, BKP-810IRB, BKPT-1216IRB, BKPL-1822IRB
Description:	SealTech Resealable Device Pouch is a compact, FDA-compliant clear pouch made from LDPE/Nylon, featuring a convenient resealable closure; irradiated. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility 10 ⁻⁶ Sterility Assurance Level in compliance with AAMI Guidelines for Sterilization.
Physical Data:	
Materials of Construction:	Clear Laminate: (LDPE) Polyethylene 75%/Cast Nylon 25%
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Thickness: 3 mil (0.076 mm) Open Width x Height: BKPH-58IRB: 5.5 in x 9 in (13.97 cm x 22.86 cm) BKP-810IRB: 8 in x 10 in (20.32 cm x 25.4 cm) BKPT-1216IRB: 12 in x 16 in (30.48 cm x 40.64 cm) BKPL-1822IRB: 18 in x 22 in (45.72 cm x 55.88 cm)
For Use with: (as needed)	BKPH-58IRB, BKP-810IRB: Mobile Phones/Smart Phones BKPT-1216IRB: Tablets/Touchscreen Devices BKPL-1822IRB: Laptops
Autoclavability:	Non-Autoclavable
Chemical Compatibility:	Good compatibility with standard disinfectants and cleaning solutions; avoid chlorites.

Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	5 pouches per bag; 10 pouches per double bagged pack; 50 pouches per case.
Traceability:	Lot number and part number printed on each package and each case.
Shelf Life:	36 months
Irradiation Standards:	Minimum Irradiation dose 25 kGy
Certificates:	Certificate of Conformance and Certificate of Processing shipped with each order and available on website. Validation paperwork available upon request.