

V-BLADE MOP, NOVAPOLY

Product Specification

Product Name:	V-Blade Mop, NovaPoly
Manufacturer Part Number(s):	VNTC6-28D, VNTC6-42D
Description:	<p>A slip-on smooth polyester mop cover to be used with the dual blade QPVNT-28 or QPVNT-42 series mop adapters. Used for cleaning and sanitizing in and around process and lab fixtures in clean and controlled environments.</p> <p>Manufactured at Micronova's Torrance, CA facility. Designed for cleaning and sanitizing in and around process and lab equipment in clean and controlled environments.</p>
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer Layer:	100% Smooth polyester
Interior:	Thin urethane foam
Dimensions:	
Folded for use:	VNTC6-42D: 4 in x 2 in (10 cm x 5 cm) VNTC6-28D: 2 in x 8 in (5 cm x 20 cm)
Opened flat:	VNTC6-42D: 4 in x 4 in (10 cm x 10 cm) VNTC6-28D: 2 in x 16 in (5 cm x 40 cm)
For Use with: (as needed)	<p>QPVNT-28 Stainless Steel and QPVNT-42 Plastic V-Blade adapters. V-Blade adapter attaches to the STA series threaded handles.</p> <p>Match cover suffix numbers with adapter suffix numbers.</p>
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	Individually bagged 24 per case
Traceability:	Part number and lot number on each bag and case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request