

SLIM-T AND BOX WRINGERS

Product Specification

Product Name:	Slim-T and Box Wringers
Manufacturer Part Numbers:	W-9, W-92, W-9B and W-9BW
Description:	Simple electropolished stainless steel mop wringer. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Clips:	Electropolished Stainless Steel with holes Nylon: W-9BW includes nylon clips to attach to B-99 buckets
Dimensions:	W-9: 18 in x 5 in x 5 in – 1 1/2 in Lip W-92: 16 in x 5 in x 4 1/2 in – 4 1/2 in Lip W-9B: 10 1/2 in x 9 1/2 in x 4 1/4 in – Cube with angled base W-9BW: 10 1/2 in x 9 1/2 in x 4 1/4 in – Cube with angled base and Nylon Clips to attach to B-99 Buckets
For Use with: (as needed)	W-9, W-92: Wring flat head or T-Mops and B-7 series or B-99 buckets. W-9B: String style mops into B-7 bucket. W-9BW: String style mops into B-99 bucket.
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions
Packaging:	Individually bagged
Traceability:	Lot Number and part number on each bag and carton
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request