

# BioGuard™ Premium Biohazard Bag, Irradiated

## Product Specification

Product Name:	BioGuard™ Premium Biohazard Bag, Irradiated
Manufacturer Part Number(s):	BHPR-21REIRBB
Description:	BioGuard™ Premium Biohazard Liner, Red, 21 X 36 inches with 70% IPA/Isopropyl Alcohol resistant black biohazard ink. Irradiated. Does not include pre-printed lines for user information.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Laminate: 28% METPET 72% LDPE
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Thickness: 3.6 mil (0.09 mm) Open Width x Height: 21 in x 36 in (534 mm x 916 mm)
For Use with: (as needed)	N/A
Autoclavability:	Not autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	2 per pack, double bagged, 25 packs of 2 each per case
Traceability:	Lot number and part number printed on each package and each case.
Shelf Life:	36 months
Irradiation Standards:	Minimum Irradiation dose 25 kGy
Certificates:	Certificate of Conformance and Certificate of Processing shipped with each order. Also available upon request and on the website.