

PLASTIC TUBE BRUSHES

Product Specification

Product Name:	Plastic Tube Brushes
Manufacturer Part Number(s):	BRPL-1418, BRPL-1218, BRPL-3418, BRPL-118, BRPL-25S36, BRPL-32548 BRPL-1218H, BRPL-160H, BRPL-236H, BRPL-237520H
Description:	Tube brushes constructed from polypropylene with microfilament fiber brush heads and flexible handles. Available with optional integrated handle holes to facilitate hanging; part numbers with H at the end. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Handle: Polypropylene Bristle: Polypropylene with microfilament fibers
Packaging Materials:	Packaged in LDPE bag, LDPE case liner.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Standard Brushes: Brush Face: 4 in (10.16 cm) Overall Length: 18 in (46 cm) Brush Diameter: BRPL-1418: ¼ in (.64 cm) BRPL-1218: ½ in (1.27 cm) BRPL-3418: ¾ in (1.91 cm) BRPL-118: 1 in (2.54 cm) Standard Brushes with Handle Holes: Brush Face: 4 in (10.16 cm) Overall Length: 18 in (46 cm) Brush Diameter: BRPL-1218H: ½ in (1.27 cm) Extended Length Brushes: Brush Face: 5 in (12.7 cm) Brush Diameter: BRPL-25S36: 2 in (5 cm) Overall Length: 36 inches (91.44 cm)

	BRPL-32548: 3 ¼ in (8.26 cm) Overall Length: 48 inches (121 cm)
	Extended Length Brushes with Handle Holes: Brush Diameter:
	BRPL-160H: 1 in (2.54 cm) Brush Face: 5 in (12.7 cm) Overall Length: 60 inches (152.4 cm)
	BRPL-236H: 2 in (5 cm) Brush Face: 4 in (10.16 cm) Overall Length: 36 inches (91.44 cm)
	BRPL-237520H: 2 ¾ in (2.03 cm) Brush Face: 5 in (12.7 cm) Overall Length: 20 inches (50.8 cm)
For Use with: (as needed)	N/A
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Storage Conditions:	Ambient temperature (Recommended between 15°C - 30°C.) Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Packaging Configuration:	Each brush is individually packaged
Traceability:	Every bag is labeled with its lot number and part number.
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on the website.