BUCKETLINER™

Product Specification

Product Name:	BucketLiner™
Manufacturer Part Number(s):	BKL7-3119, BKL-2424, BKL1-1818, BKL99-3521
Description:	Low-density polyethylene/nylon laminate liners are made from FDA-compliant LDPE, with a flat seal bottom, designed to provide an impermeable barrier for lining buckets and containers. The liners can be used for liquids or dry materials. Manufactured at Micronova's facilities in USA. Designed for cleanroom and controlled environments. ISO 5 to ISO 9
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Clear Laminate: FDA-compliant LDPE Polyethylene 75%/Cast Nylon 25%
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Thickness: 3 mil (0.076 mm) Open Width x Height BKL1-1818: 18 x 18 in (45.7 x 45.7 cm) BKL-2424: 24 x 24 in (61 x 61 cm) BKL7-3119: 31 x 19 in (78.7 x 48.3 cm) BKL99-3521: 35 x 21 in (89 x 53.3 cm)

Micronova Manufacturing
3431 Lomita Boulevard
Torrance, CA 90505
Phone: 310-784-6990 / micronova-mfg.com

For Use with: (as needed)	Buckets and containers similar in shape to: BKL1-1818: 1-gal cylindrical or square pails BKL-2424: 6-gal cylindrical or square bucket BKL7-3119: 6-gal rectangular bucket; B-7 Slim T™ series buckets BKL99-3521: 8-gal rectangular bucket; B-99 stainless steel bucket
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions; avoid chlorites.
Storage Conditions:	Ambient temperature (recommended 15°C - 30°C or 59°F - 86°F). Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet light sources.
Packaging Configuration:	BKL1-1818: 50 liners per pack; 3 packs per case; 150 liners total. BKL7-3119, BKL-2424, BKL99-3521: 20 liners per pack; 10 packs per case; 200 liners total.
Traceability:	Every bag and case is labeled with its part number and lot number.
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on the website.

Micronova Manufacturing
3431 Lomita Boulevard
Torrance, CA 90505
Phone: 310-784-6990 / micronova-mfg.com