AUTOCLAVE INDICATOR TAPE, LEAD FREE

Product Specification

Product Name:	Autoclave Indicator Tape, Lead Free
Manufacturer Part Number(s):	AUTOPC-34, AUTOPC-12, AUTOPC-1, AUTOPC-15, AUTOPC-2, AUTOPC-25, AUTOPC-3, AUTOPC-1PK, AUTOPC-15PK, AUTOPC-2PK, AUTOPC-3PK
Description:	A lead-free, latex-free, indicator tape for steam sterilization cycles. Green indicator lines change to black at completion of the sterilization cycle in both gravity discharge and vacuum autoclave. Available in a compact 6-roll pack size, identified by product numbers ending in "PK." Manufactured at Micronova's facilities in USA. Designed for cleanroom and controlled environments. ISO 5 to ISO 9
Manufacturing Standards:	ISO 9001 Certified Facility RoHS Compliant International Standard per ISO 11140-1:2014
Physical Data:	
Materials of Construction:	Film: Saturated crepe paper Backing: Synthetic, solvent-free adhesive Core: Plastic Color: Natural with green indicator lines
Packaging Materials:	Packaged in LDPE bag, LDPE case liner. ECT rated case of recycled material.
BSE/TSE Statement:	We certify that this product is manufactured without the use of raw materials of animal or human origin.
Dimensions:	Core: 3 inch (7.6 cm) Length: 60 yards (55 m) Width: AUTOPC-34: ¾ in (1.9 cm) AUTOPC-12: ½ in (1.27 cm) AUTOPC-1: 1 in (2.5 cm) AUTOPC-15: 1.5 in (3.8 cm) AUTOPC-2: 2 in (5.7 cm) AUTOPC-25: 2.5 in (6.4 cm) AUTOPC-3: 3 in (7.6 cm)

Micronova Manufacturing
3431 Lomita Boulevard
Torrance, CA 90505
Phone: 310-784-6990 / micronova-mfg.com

	AUTOPC-1PK: 1 in (2.5 cm) AUTOPC-15PK: 1.5" (3.8 cm) AUTOPC-2PK: 2 in (5.7 cm) AUTOPC-3PK: 3 in (7.6 cm)
For Use with: (as needed)	N/A
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with standard disinfectants and cleaning solutions.
Storage Conditions:	Ambient temperature (recommended 15°C - 30°C or 59°F - 86°F) in its original packaging/container. Avoid direct sunlight and store in an area protected from extreme temperature, humidity, and UV exposure. Keep away from chemicals and cleaning agents. Store on a flat surface to prevent tape telescoping.
Packaging Configuration:	Individually bagged- AUTOPC-1: 36 rolls per case AUTOPC-12; AUTOPC-34: 16 rolls per case AUTOPC-15; AUTOPC-2: 24 rolls per case AUTOPC-25; AUTOPC-3: 12 rolls per case Small pack; Individually bagged- AUTOPC-1PK, AUTOPC-15PK, AUTOPC-2PK, AUTOPC-3PK: 6 rolls per pack
Traceability:	Every bag and case is labeled with its lot number and part number.
Shelf Life:	N/A
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website.

Micronova Manufacturing
3431 Lomita Boulevard
Torrance, CA 90505
Phone: 310-784-6990 / micronova-mfg.com